

PHILIPS

EXPERIENCE GUIDANCE LIKE NEVER BEFORE

TRACKING BRACKETS



Accurate Guidance. Confident Outcomes.

Tracking Technology

Tracking brackets are designed with a sensor mount to hold customized electromagnetic tracking sensors (EMT) in close proximity to the ultrasound transducer while connected to your ultrasound system.

- When ready for puncture, a disposable needle guide is used, providing a path for directing instruments according to on-screen software guidelines.
- Included with tracking bracket are five sterile procedure kits including: needle guide, transducer cover, gel packet, elastic bands and non-sterile probe clips
- General Purpose Electromagnetic Sensors are available by contacting your Philips representative.

PercuNav Tracking Bracket Starter Kits

C5-1



Ultra-Pro II™ style bracket
REF. #667-130
Replacement Kit 610-608
Lock Pin Kit 610-1007

X6-1



Ultra-Pro II™ style bracket
REF. #667-131
Replacement Kit 610-608
Lock Pin Kit 610-1007

C9-2



Ultra-Pro II™ style bracket
REF. #667-138
Replacement Kit 610-608
Lock Pin Kit 610-1007

C10-4ec



REF. #667-147
Replacement Kit 667-134

L12-5 50mm



Infiniti Plus™ style bracket
REF. #667-132
Replacement Kit 610-1232

*To order General Purpose Electromagnetic Sensors contact Philips

*To order Tracking Bracket contact CIVCO: 800-445-6741

Notes:





Accurate Guidance. Confident Outcomes.

For international distribution, visit CIVCO Worldwide at www.CIVCO.com

Worldwide Headquarters:

102 First Street South, Kalona, Iowa 52247-9589 USA

319/656-4447, 800/445-6741 • Fax: 319/656-4451, 877/329-2482

www.CIVCO.com

COPYRIGHT © 2014. CIVCO IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE THE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. PRINTED IN USA. 2014R-2809