

The Right Way

to Store Probes & Stay Compliant

“Store the device in a manner that will protect from damage or contamination... Hanging vertically in a cabinet and storing in a clean environment.”
 -The Joint Commission



Safe & Clean Storage For Disinfected Ultrasound Probes

CIVCO's Storage Systems provide a safe and clean storage environment for your disinfected transesophageal and general purpose/endocavity ultrasound probes.

- Fan draws air through replaceable HEPA filter for continuous clean filtered air flow
- Filter life is monitored by the storage system and an indicator light will signal when filter replacement is required
- Lockable door for secure and safe storage

REF #	PRODUCT DESCRIPTION	QTY
610-1323	TEE Probe Cabinet; Holds up to 3 probes	1
610-1327	General Purpose & Endocavity Probe Cabinet; Holds up to 6 probes	1
610-1328	Storage Cabinet replacement HEPA filter	1



CIVCO.com | 800.445.6741 | 319.248.6757